APCapplicator

DE	Argonplasma-Applikator	3
EN	Argon Plasma Applicator	21
USA	Argon Plasma Applicator	37
FR	Applicateur de plasma d'argon	55
ES	Aplicador de plasma de argón	73
IT	Manipolo a plasma d'argon	91
PT	Aplicador de plasma de argônio	109
EL	Εφαρμογέας πλάσματος αργού	127
NL	Argonplasma-applicator	145
DA	Argonplasma-applikator	163
SV	Argonplasma-applikator	181
FI	Argonplasma-applikaattori	197
PL	Aplikator plazmy argonowej	213
CS	Aplikátor argonové plazmy	231
HU	Argonplazma-applikátor	247
RU	Аргоноплазменный аппликатор	263
TR	Argon plazma aplikatörü	281
ZH	氩气电极	297
KO	아르곤 플라즈마 어플리케이터	311

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DE

APCapplicator

Argonplasma-Applikator

20132-250, 20132-251, 20132-252, 20132-253, 20132-254,
20132-255, 20132-256, 20132-257

STERILE EO

Inhalt

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
STERILE EO

Table of Contents

1	Product description	23
2	General instructions for use	23
2.1	Intended use / Indications for use	24
2.2	Compatibility	24
2.3	Intended patient population	25
2.4	Contraindications	25
2.5	Side effects	25
2.6	Single-use	25
2.7	Environment	25
2.8	Qualification of user	25
2.9	Performance characteristics	26
2.10	Overall length	26
3	Maximum electrical capacity	26
4	Safety instructions	26
5	Required accessories	29
6	How to use	30
6.1	Checking the expiration date	30
6.2	Checking and opening the packaging	30
6.3	Checking the product	30
6.4	Connecting the product and selecting unit settings	30
6.5	In the case of endoscopic application: Selecting a Trocár	33
6.6	Using the product	34
7	Disposal	35
8	Symbols	35

1 Product description

- 1 Yellow button: activates CUT
- 2 Blue button: activates COAG
- 3 Light gray ReMode button: Toggles between two setting options of a program
- 4 Rotary wheel: For freely adjustable movement of the HF electrode
- 5 Shaft tube: For turning spatula electrodes
- 6 APC connector: For connection to the APC unit

2 General instructions for use

These notes on use describe the normal use of the product. Please read all information carefully.

These notes on use do not replace the user manuals of the units and the products that are used with the instrument! Read the user manuals of the units and products and ask Erbe or your distributor in case of doubt!

Note: Report serious incidents with the product to your local dealer or Erbe. If you are a user in the European Union, also report incidents to the responsible authority in your Member State.

2.1 Intended use / Indications for use

20132-256, 20132-257: The Argon Plasma Applicators are intended for monopolar coagulation and devitalization of tissue using argon plasma. They are used in open surgery.

20132-250, 20132-251, 20132-252, 20132-253, 20132-254 and 20132-255: The Argon Plasma Applicators are intended for monopolar cutting and coagulation of tissue with and without argon gas or for coagulation and devitalization using argon plasma:

- used in open surgery 20132-250, 20132-251, 20132-252, 20132-253
- used in laparoscopic surgery 20132-254, 20132-255.

2.2 Compatibility

The Argon Plasma Applicator is intended for use with the following units:

- APC 2 in combination with an electrosurgical unit:
 - VIO 200 D, VIO 300 D
 - VIO 200 S, VIO 300 S
- APC 3 in combination with a VIO 3 or VIO 3n electrosurgical unit

The adapter REF No. 20132-249 is required for connection to an APC socket with separate connections for HF and argon gas.

2.3 Intended patient population

No restrictions.

2.4 Contraindications

No known contraindications.

2.5 Side effects

Known side effects are:

- neuromuscular stimulation
-

2.6 Single-use

The product is intended for single use.

The product is supplied sterile.

2.7 Environment

For the intended use, the product may only be used in premises for medical purposes.

2.8 Qualification of user

For the intended use, the product should only be used by medical professionals who have been trained on the basis of these notes on use.

2.9 Performance characteristics

Performance characteristics relating to the intended use are:

- Conducting high-frequency current from the connector to the distal end.
- Passing argon gas from the connector to the distal end.

2.10 Overall length

The overall length of the product (including connector, cable/tube) is 3.7 m.

3 Maximum electrical capacity

The maximum electrical capacity of this product is:

- 5000 V_p

4 Safety instructions

WARNING!

The activated electrosurgical unit can impair the performance of or damage cardiac pacemakers, internal defibrillators, and other active implants. If patients have active implants, consult the manufacturer of the implant or the relevant medical specialist prior to procedure.

The activated instrument may damage or influence the functionality of passive implants or damage metallic objects in or on the patient's body. Be careful when activating in the vicinity of passive implants or metallic objects.

Check the product for visible damage before use. Do not use this product if it is damaged.

Do not use in the presence of combustible or explosive materials!

Protect this product from any mechanical damage! Do not throw! Do not use force!

Do not kink electric connecting cables and/or tubing, only wrap them up loosely if necessary. Do not wind around objects.

Inadvertently activated instruments and instruments that are still hot after use can cause burns or fires:

- Store instrument away from the patient, never on the patient.
- Avoid contacting hot instrument components (e.g. the distal end).
- Never set the instrument down on flammable materials.

Risk of burns. Do not activate the instrument while cleaning.

Nerves may become thermally damaged or electrically stimulated when the product is used in close proximity to nerve structures. This may cause spasms and muscle contractions. Risk of injury to the patient. Use the lowest possible unit settings to achieve the desired tissue effects.

Flammable endogenous gases in the gastrointestinal tract. Risk of explosion for the patient! Risk of burns for the user. Follow therapeutic guidelines for bowel preparation. Mitigate the presence of such gases before using the instrument and/or insufflate with CO₂.

To minimize gas emphysema and gas embolisms, use non-contact technique and the lowest possible argon flow rate to achieve the desired tissue effect. Be aware of open vessels during APC activation.

Argon gas may not immediately escape from the anatomical application area. Be aware of signs of increased intraluminal gas pressure, for example, by palpation and/or visual inspection. Extract gas as needed.

Use of the instrument under too high oxygen concentration in the tracheobronchial system. Risk of fires in the tracheobronchial system for the patient! Reduce oxygen concentration according to the therapeutic guidelines before activating the instrument in the tracheobronchial system.

Avoid activation on conductive materials, e.g. metal, as this may cause unintended thermal insult.

Risk of burns as a result of leakage current. Position the electric connecting cables of the instrument so that they do not touch the patient and are as far apart as possible from other electric cables.

Risk of capacitive coupling. Electrosurgical current can be transferred from the activated instrument to other conductive materials even if the two are not touching ("capacitive coupling"). Risk of burns to the patient and/or medical personnel, e.g.:

- Position the electric connecting cables of the instruments as far apart as possible.
- Place the instrument at the usual working distance from the tissue before activating the instrument.
- Do not use hybrid trocars (e.g., a combination of metal and plastic). Only use trocars that can be used with electrosurgical accessories.
- Only use endoscopes that can be used for electrosurgery. Please consult the endoscope manufacturer.

If the instrument is activated for a long time and/or activated in quick succession (with no cooling time), the distal instrument tip can become extremely hot. This is not a fault in the material but rather a consequence of heavy use. Excessive heat on the distal end of the instrument can also cause unwanted tissue damage in the patient; therefore, avoid activating the instrument for long periods.

Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in surgical applications, e.g. in open and laparoscopic procedures.

The product must not be processed. Processing can affect the material properties and/or function of the product such that application in accordance with the intended use is no longer possible. It also cannot be ensured that the product will be sterile after processing. Risks include:

- Components/connections that are damaged by high temperatures (e.g., in the autoclave)
- Components (e.g., filter, extremely small lumen, etc.) that may prevent internal cleaning/disinfection

Erbe Elektromedizin expressly warns against modifying the product. Any modification exempts Erbe Elektromedizin from any and all liability.

5 Required accessories

Note: The APC units are equipped with different APC sockets. Depending on the type of socket, you may or may not need an adapter to connect the applicator.

Socket type	Accessories for connecting the applicator
	No other accessories are required.
	Adapter REF No. 20132-249

6 How to use

6.1 Checking the expiration date

The expiration date is indicated on the packaging. Do not use this product if the expiration date has passed!

6.2 Checking and opening the packaging

WARNING! Only use the product if the sterile packaging is intact and unopened. Damaged packaging, and potentially non sterile packaging can be identified for example by damaged seal seams, illegible labeling, damp or even wet areas, etc.

Use aseptic techniques when removing the product from the packaging!

6.3 Checking the product

WARNING! Check the product for visible damage before use. Do not use this product if it is damaged.

6.4 Connecting the product and selecting unit settings

WARNING! Risk of injury! Do not activate the electrosurgical unit while connecting the instrument.

WARNING! Select unit settings that are as low as possible to achieve the desired tissue effect.


Only select unit settings that do not exceed the maximum electrical capacity of the product (V_p max or kV_p max) as this could damage the product. If several products (e.g., adapter and probe) with different electrical capacities are combined, do not exceed the lowest electrical capacity. For information on how

to configure suitable unit settings, please refer to the user manual of the Erbe electrosurgical unit (see "Accessories" section).

Note: The APC units are equipped with different APC sockets. Depending on the type of socket, you may or may not need an adapter to connect the applicator.


To determine your type of socket please read the chapter "Required accessories".

Connection to an APC 2 with APC socket


Note: This type of socket is available in two variants, with and without finger switch activation. For the applicator, you will need the variant with finger switch activation, which can be recognized by the large metal ring at the left contact of the socket.

1. Switch on the VIO electrosurgical unit and select the basic program.
2. Insert the APC connector (1) into the APC socket (2).


3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 2 with APC socket

Note: The adapter REF No. 20132-249 must be plugged into the APC socket!

1. Switch on the VIO electrosurgical unit and select the basic program.
2. Insert the APC connector (1) into the adapter (2).

3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 3:

1. Switch on the VIO electrosurgical unit and pull up the main menu.
2. Insert the APC connector into the APC socket on the APC 3.

3. For the connected applicator, a standard setting for mode, effect level and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

6.5 In the case of endoscopic application: Selecting a Trocar

Select a trocar for instruments with diameter 5 mm.

Flash-overs, activation in the trocar or without visual control of the distal end of the instrument can cause unwanted thermal injuries to the patient. Therefore, check that the distal end of the instrument protrudes at least 10 mm from the trocar before each activation.

WARNING! The HF leakage current at the ocular of the endoscope may cause thermal injuries to the user. Only use endoscopes that can be used for electrosurgery. Please consult the endoscope manufacturer.

WARNING! Unwanted stimulation of nerves and muscles due to low-frequency currents. Dangerous summation of patient leakage currents from several power sources: It is not possible to entirely avoid low-frequency currents in electrosurgery, but they can be limited by applying technical measures to equipment and accessories. Limit types show how wide-ranging these measures are, e.g., type B, BF, or CF as per EN 60601. Type CF offers the highest level of protection against low-frequency currents. However, please note the following: If you use an electrosurgical unit and an endoscope that is supplied with energy at the same time, the low-frequency currents from both power sources accumulate to form one single patient leakage current. If the electrosurgical unit belongs to type CF and the endoscope to type BF, the less favorable limit values for type BF apply to the overall patient leakage current. You should therefore only use power sources with the same type classification, ideally those of type CF.

6.6 Using the product

WARNING! Do not bend the shaft tube in the case of applicators with a rotary wheel and movable electrode.

WARNING! The applicator with a malleable shaft (REF No. 20132-256, 20132-257) should only be bent manually, do not use any tools.

CAUTION! Biological material can impair the effectiveness of the instrument. Keep the instrument clean during use. Do not activate the instrument while cleaning.

Use a non-abrasive material, e.g. gauze and sterile water as necessary for cleaning.

In the case of endoscopic application: Retract the electrode completely when introducing the applicator or removing it from the trocar. You must feel the rotary wheel clicking into place. Only introduce the applicator into the trocar under visual control.

Press the yellow button to activate cutting mode (CUT). Press the blue button to activate coagulation mode (COAG). Make sure that the electrosurgical unit responds correctly.

Using the light gray ReMode button, you can, in the case of a program with two setting options *a* and *b*, toggle between both setting options.

REF No. 20132-256, 20132-257: These applicators have only one blue button. The blue button activates the coagulation mode. Make sure that the electrosurgical unit responds correctly.

In order to move the electrode longitudinally, turn the rotary wheel, or alternatively, move the shaft tube back and forth.

Spatula electrodes can be turned either way in any position. To do so, turn the shaft tube.

WARNING! Instrument tips that are pressed to tissue could cause a mechanical perforation. Be careful with contact between the instrument tip and tissue and avoid unintentional contact with tissue.

7 Disposal

WARNING! Risk of infection by contaminated and potentially infectious material, e.g. used instrument, etc.

Dispose of the product, packaging material, and accessories (if available) in accordance with the applicable guidelines and laws.

8 Symbols

Symbol	Explanation	Symbol	Explanation
	Consult instructions for use		Caution, read safety instructions before use.
	Catalogue number		Batch code
	Manufacturer		Keep away from sunlight
	Keep dry		Medical device
	Quantity (x)		Use-by date ¹⁾

Symbol	Explanation	Symbol	Explanation
	Sterilized using ethylene oxide		Do not use if package is damaged and consult user manual
	Do not reuse		Rigid shaft tube with needle
	Rigid shaft tube with scoop		Malleable shaft tube
	Single sterile barrier system ²⁾		Do not resterilize
	European conformity marking		

- 1.) The date next to the symbol is in the format Year-Month-Day.
- 2.) The date next to the symbol states the Year-Month-Day on which the sterile barrier system was sealed.

NOTES ON USE

USA

APCapplicator

Argon Plasma Applicator

20132-250, 20132-251, 20132-252, 20132-253, 20132-254,
20132-255, 20132-256, 20132-257


STERILE EO

Table of Contents

1	Product description	39
2	General instructions for use	39
2.1	Intended use / Indications for use	40
2.2	Compatibility	40
2.3	Intended patient population	41
2.4	Contraindications	41
2.5	Side effects	41
2.6	Single-use	41
2.7	Environment	41
2.8	Qualification of user	41
3	Maximum electrical capacity	41
4	Safety instructions	42
5	Required accessories.	45
6	How to use	46
6.1	Checking the expiration date.	46
6.2	Checking and opening the packaging	46
6.3	Checking the product.	46
6.4	Connecting the product and selecting unit settings	46
6.5	In the case of endoscopic application: Selecting a Trocár.	49
6.6	Using the product.	50
7	Disposal	51
8	Symbols	51

1 Product description

- 1 Yellow button: activates CUT
- 2 Blue button: activates COAG
- 3 Light gray ReMode button: Toggles between two setting options of a program
- 4 Rotary wheel: For freely adjustable movement of the HF electrode
- 5 Shaft tube: For turning spatula electrodes
- 6 APC connector: For connection to the APC unit

2 General instructions for use

These notes on use describe the normal use of the product. Please read all information carefully.

These notes on use do not replace the user manuals of the units and the products that are used with the instrument! Read the user manuals of the units and products and ask Erbe or your distributor in case of doubt!

CAUTION! Federal law (USA) restricts this device to sale by or on the order of a physician!

Please observe the warnings in chapter "Safety instructions".

2.1 Intended use / Indications for use

20132-256, 20132-257: The Argon Plasma Applicators are intended for monopolar coagulation and ablation of tissue using argon plasma. They are used in open surgery.

20132-250, 20132-251, 20132-252, 20132-253, 20132-254 and 20132-255: The Argon Plasma Applicators are intended for monopolar cutting and coagulation of tissue with and without argon gas or for coagulation and ablation using argon plasma:

- used in open surgery 20132-250, 20132-251, 20132-252, 20132-253
- used in laparoscopic surgery 20132-254, 20132-255

2.2 Compatibility

The APCApplicator is intended for use with the following units:

- APC 2 in combination with the following electrosurgical units:
 - VIO 300 D
 - VIO 200 S, VIO 300 S
- APC 3 in combination with a electrosurgical unit VIO 3

The adapter 20132-249 is required for connection to an APC socket with separate connections for HF and argon gas.

2.3 Intended patient population

No restrictions.

2.4 Contraindications

No known contraindications.

2.5 Side effects

Known side effects are:

- neuromuscular stimulation

2.6 Single-use

The product is intended for single use.

The product is supplied sterile.

2.7 Environment

For the intended use, the product may only be used in premises for medical purposes.

2.8 Qualification of user

For the intended use, the product should only be used by medical professionals who have been trained on the basis of these notes on use.

3 Maximum electrical capacity

The maximum electrical capacity of this product is:

- 5000 V_p

4 Safety instructions

WARNING!

The activated electrosurgical unit can impair the performance of or damage cardiac pacemakers, internal defibrillators, and other active implants. If patients have active implants, consult the manufacturer of the implant or the relevant medical specialist prior to procedure.

The activated instrument may damage or influence the functionality of passive implants or damage metallic objects in or on the patient's body. Be careful when activating in the vicinity of passive implants or metallic objects.

Check the product for visible damage before use. Do not use this product if it is damaged.

Do not use in the presence of combustible or explosive materials!

Protect this product from any mechanical damage! Do not throw! Do not use force!

Do not kink electric connecting cables and/or tubing, only wrap them up loosely if necessary. Do not wind around objects.

Inadvertently activated instruments and instruments that are still hot after use can cause burns or fires:

- Store instrument away from the patient, never on the patient.
- Avoid contacting hot instrument components (e.g. the distal end).
- Never set the instrument down on flammable materials.

Risk of burns. Do not activate the instrument while cleaning.

Nerves may become thermally damaged or electrically stimulated when the product is used in close proximity to nerve structures. This may cause spasms and muscle contractions. Risk of injury to the patient. Use the lowest possible unit settings to achieve the desired tissue effects.

Flammable endogenous gases in the gastrointestinal tract. Risk of explosion for the patient! Risk of burns for the user. Follow therapeutic guidelines for bowel preparation. Mitigate the presence of such gases before using the instrument and/or insufflate with CO₂.

To minimize gas emphysema and gas embolisms, use non-contact technique and the lowest possible argon flow rate to achieve the desired tissue effect. Be aware of open vessels during APC activation.

Argon gas may not immediately escape from the anatomical application area. Be aware of signs of increased intraluminal gas pressure, for example, by palpation and/or visual inspection. Extract gas as needed.

Use of the instrument under too high oxygen concentration in the tracheobronchial system. Risk of fires in the tracheobronchial system for the patient! Reduce oxygen concentration according to the therapeutic guidelines before activating the instrument in the tracheobronchial system.

Avoid activation on conductive materials, e.g. metal, as this may cause unintended thermal insult.

Risk of burns as a result of leakage current. Position the electric connecting cables of the instrument so that they do not touch the patient and are as far apart as possible from other electric cables.

Risk of capacitive coupling. Electrosurgical current can be transferred from the activated instrument to other conductive materials even if the two are not touching ("capacitive coupling"). Risk of burns to the patient and/or medical personnel, e.g.:

- Position the electric connecting cables of the instruments as far apart as possible.
- Place the instrument at the usual working distance from the tissue before activating the instrument.
- Do not use hybrid trocars (e.g., a combination of metal and plastic). Only use trocars that can be used with electrosurgical accessories.
- Only use endoscopes that can be used for electrosurgery. Please consult the endoscope manufacturer.

If the instrument is activated for a long time and/or activated in quick succession (with no cooling time), the distal instrument tip can become extremely hot. This is not a fault in the material but rather a consequence of heavy use. Excessive heat on the distal end of the instrument can also cause unwanted tissue damage in the patient; therefore, avoid activating the instrument for long periods.

Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in surgical applications, e.g. in open and laparoscopic procedures.

The product must not be processed. Processing can affect the material properties and/or function of the product such that application in accordance with the intended use is no longer possible. It also cannot be ensured that the product will be sterile after processing. Risks include:

- Components/connections that are damaged by high temperatures (e.g., in the autoclave)
- Components (e.g., filter, extremely small lumen, etc.) that may prevent internal cleaning/disinfection

Erbe Elektromedizin expressly warns against modifying the product. Any modification exempts Erbe Elektromedizin from any and all liability.

5 Required accessories

Note: The APC units are equipped with different APC sockets. Depending on the type of socket, you may or may not need an adapter to connect the applicator.

Socket type	Accessories for connecting the applicator
	No other accessories are required.
	Adapter REF No. 20132-249

6 How to use

6.1 Checking the expiration date

The expiration date is indicated on the packaging. Do not use this product if the expiration date has passed!

6.2 Checking and opening the packaging

WARNING! Only use the product if the sterile packaging is intact and unopened. Damaged packaging, and potentially non sterile packaging can be identified for example by damaged seal seams, illegible labeling, damp or even wet areas, etc.

Use aseptic techniques when removing the product from the packaging!

6.3 Checking the product

WARNING! Check the product for visible damage before use. Do not use this product if it is damaged.

6.4 Connecting the product and selecting unit settings

WARNING! Risk of injury! Do not activate the electrosurgical unit while connecting the instrument.

WARNING! Select unit settings that are as low as possible to achieve the desired tissue effect.


Only select unit settings that do not exceed the maximum electrical capacity of the product (V_p max or kV_p max) as this could damage the product. If several products (e.g., adapter and probe) with different electrical capacities are combined, do not exceed the lowest electrical capacity. For information on how

to configure suitable unit settings, please refer to the user manual of the Erbe electrosurgical unit (see "Accessories" section).

Note: The APC units are equipped with different APC sockets. Depending on the type of socket, you may or may not need an adapter to connect the applicator.


To determine your type of socket please read the chapter "Required accessories".

Connection to an APC 2 with APC socket


Note: This type of socket is available in two variants, with and without finger switch activation. For the applicator, you will need the variant with finger switch activation, which can be recognized by the large metal ring at the left contact of the socket.

1. Switch on the VIO electrosurgical unit and select the basic program.
2. Insert the APC connector (1) into the APC socket (2).


3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 2 with APC socket

Note: The adapter REF No. 20132-249 must be plugged into the APC socket!

1. Switch on the VIO electrosurgical unit and select the basic program.
2. Insert the APC connector (1) into the adapter (2).

3. For the connected applicator, a standard setting for mode, power, effect, and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

Connection to an APC 3:

1. Switch on the VIO electrosurgical unit and pull up the main menu.
2. Insert the APC connector into the APC socket on the APC 3.

3. For the connected applicator, a standard setting for mode, effect level and gas flow is automatically made on the VIO electrosurgical unit. You can change this standard setting within a certain range as required.

6.5 In the case of endoscopic application: Selecting a Trocar

Select a trocar for instruments with diameter 5 mm.

Flash-overs, activation in the trocar or without visual control of the distal end of the instrument can cause unwanted thermal injuries to the patient. Therefore, check that the distal end of the instrument protrudes at least 10 mm from the trocar before each activation.

WARNING! The HF leakage current at the ocular of the endoscope may cause thermal injuries to the user. Only use endoscopes that can be used for electrosurgery. Please consult the endoscope manufacturer.

WARNING! Unwanted stimulation of nerves and muscles due to low-frequency currents. Dangerous summation of patient leakage currents from several power sources: It is not possible to entirely avoid low-frequency currents in electrosurgery, but they can be limited by applying technical measures to equipment and accessories. Limit types show how wide-ranging these measures are, e.g., type B, BF, or CF as per EN 60601. Type CF offers the highest level of protection against low-frequency currents. However, please note the following: If you use an electrosurgical unit and an endoscope that is supplied with energy at the same time, the low-frequency currents from both power sources accumulate to form one single patient leakage current. If the electrosurgical unit belongs to type CF and the endoscope to type BF, the less favorable limit values for type BF apply to the overall patient leakage current. You should therefore only use power sources with the same type classification, ideally those of type CF.

6.6 Using the product

WARNING! Do not bend the shaft tube in the case of applicators with a rotary wheel and movable electrode.

WARNING! The applicator with a malleable shaft (REF No. 20132-256, 20132-257) should only be bent manually, do not use any tools.

CAUTION! Biological material can impair the effectiveness of the instrument. Keep the instrument clean during use. Do not activate the instrument while cleaning.

Use a non-abrasive material, e.g. gauze and sterile water as necessary for cleaning.

In the case of endoscopic application: Retract the electrode completely when introducing the applicator or removing it from the trocar. You must feel the rotary wheel clicking into place. Only introduce the applicator into the trocar under visual control.

Press the yellow button to activate cutting mode (CUT). Press the blue button to activate coagulation mode (COAG). Make sure that the electrosurgical unit responds correctly.

Using the light gray ReMode button, you can, in the case of a program with two setting options *a* and *b*, toggle between both setting options.

REF No. 20132-256, 20132-257: These applicators have only one blue button. The blue button activates the coagulation mode. Make sure that the electrosurgical unit responds correctly.

In order to move the electrode longitudinally, turn the rotary wheel, or alternatively, move the shaft tube back and forth.

Spatula electrodes can be turned either way in any position. To do so, turn the shaft tube.

WARNING! Instrument tips that are pressed to tissue could cause a mechanical perforation. Be careful with contact between the instrument tip and tissue and avoid unintentional contact with tissue.

7 Disposal

WARNING! Risk of infection by contaminated and potentially infectious material, e.g. used instrument, etc.

Dispose of the product, packaging material, and accessories (if available) in accordance with the applicable guidelines and laws.

8 Symbols

ISO 15223-1 Medical devices – Symbols to be used with information to be supplied by the manufacturer

Ref. no.	Symbol	Title	Explanation
5.1.1		Manufacturer	Indicates the medical device manufacturer.
5.1.4		Use-by date ¹⁾	Indicates the date after which the medical device is not to be used.
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.

Ref. no.	Symbol	Title	Explanation
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
5.2.3		Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
5.2.6		Do not resterilize	Indicates a medical device that is not to be resterilized.
5.2.8		Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
5.2.11		Single sterile barrier system	Indicates a single sterile barrier system. ²⁾
5.3.2		Keep away from sunlight	Indicates a medical device that needs protection from light sources.

Ref. no.	Symbol	Title	Explanation
5.3.4		Keep dry	Indicates a medical device that needs to be protected from moisture.
5.4.2		Do not re-use	Indicates a medical device that is intended for one single use only.
5.4.4		Caution	Read safety instructions before use.
5.7.7		Medical device	Indicates the item is a medical device.

- 1.) The date next to the symbol is in the format Year-Month-Day.
- 2.) The date next to the symbol states the Year-Month-Day on which the sterile barrier system was sealed.

Regulation (EU) 2017/745, Annex V

Symbol	Title	Explanation
	European conformity marking	Confirmation from the manufacturer that the product meets the requirements of the applicable European guidelines.

Symbols not relating to standards

Symbol	Title	Explanation
	Quantity (x)	Quantity of the product. x pieces of the product are supplied in one packaging unit.
	Rigid shaft tube with needle	This instrument has a rigid, non-bendable instrument tip with a needle electrode.
	Rigid shaft tube with scoop	This instrument has a rigid, non-bendable instrument tip with a spatula electrode.
	Malleable shaft tube	This instrument has a bendable instrument tip.

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