
Basiswissen der Plasmachirurgie

Argonplasma-Koagulation

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- 2 Eickhoff A, Repici A, Manner H, Enderle, MD. *Electrosurgical Pocket Guide for GI Interventions*; Erbe Elektromedizin GmbH
- 3 Eickhoff, A et al. *Prospective nonrandomized comparison of two modes of argon beamer (APC) tumor desobstruction: effectiveness of new pulsed APC versus forced APC*; *Endoscopy* 2007
- 4 Zenker, M. *Argon plasma Coagulation*. *GMS Krankenhhyg Interdiszip.* 2008
- 5 Taghavi SA, Soleimani SM, Hosseini-Asl SM et al. *Adrenaline injection plus argon plasma coagulation versus adrenaline injection plus hemoclips for treating high-risk bleeding peptic ulcers: a prospective, randomized trial*. *Can J Gastroenterol* 2009; 23(10): 699 – 704.
- 6 Wang HM, Hsu PI, Lo GH et al. *Comparison of hemostatic efficacy for argon plasma coagulation and distilled water injection in treating high-risk bleeding ulcers*. *J Clin Gastroenterol* 2009; 43(10): 941 – 945.
- 7 Herrera S, Bordas JM, Llach J et al. *The beneficial effects of argon plasma coagulation in the management of different types of gastric vascular ectasia lesions in patients admitted for GI hemorrhage*. *Gastrointest Endosc* 2008; 68(3): 440 – 446.

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
- 5 Taghavi SA, Soleimani SM, Hosseini-Asl SM et al. Adrenaline injection plus argon plasma coagulation versus adrenaline injection plus hemoclips for treating high-risk bleeding peptic ulcers: a prospective, randomized trial. *Can J Gastroenterol* 2009; 23(10): 699 – 704.
- 6 Wang HM, Hsu PI, Lo GH et al. Comparison of hemostatic efficacy for argon plasma coagulation and distilled water injection in treating high-risk bleeding ulcers. *J Clin Gastroenterol* 2009; 43(10): 941 – 945.
- 7 Herrera S, Bordas JM, Llach J et al. The beneficial effects of argon plasma coagulation in the management of different types of gastric vascular ectasia lesions in patients admitted for GI hemorrhage. *Gastrointest Endosc* 2008; 68(3): 440 – 446.
- 8 Manner H, Enderle MD, Pech O et al. Second-generation argon plasma coagulation: two-center experience with 600 patients. *J Gastroenterol Hepatol* 2008; 23(6): 872 – 878.
- 9 Katsinelos P, Gkagkalis S, Paroutoglou G, Chatzimavroudis G, Fasoulas K, Zavos C, Varitimadidis K, Lazaraki G, Kotronis G, Kountouras J. A prospective comparative study of blended and pure coagulation current in endoscopic mucosal resection of large sessile colorectal polyps. *Surg Laparosc Endosc Percutan Tech*. 2014 Jun; 24(3): 226 – 31.
- 10 Manner H, Rabenstein T, Pech O, Braun K, May A, Pohl J, Angelika Behrens A, Vieth M, Ell C. Ablation of residual Barrett's epithelium after endoscopic resection: a randomized long-term follow-up study of argon plasma coagulation vs. surveillance (APE study). *Endoscopy* 2014; 46(01): 6–12
- 15 Wang H, Tao M, Zhang N, Luo L, Li D, Zou H, Zhou Y, Liang S. Bronchoscopic interventions combined with percutaneous modalities for the treatment of thyroid cancers with airway invasion. *Eur Arch Otorhinolaryngol*. 2015 Feb; 272(2): 445 – 51.
- 16 Wang JW, Huang M, Zha WJ, Zhou LF, Qi X, Wang H. Flexible bronchoscopic intervention for endobronchial hamartoma. *Zhonghua Jie He He Hu Xi Za Zhi*. 2013 Dec; 36(12): 963 – 7. Chinese
- 17 Seaman JC, Musani AI. Endobronchial ablative therapies. *Clin Chest Med*. 2013 Sep; 34(3): 417 – 25. doi: 10.1016/j.ccm.2013.04.006.
- 18 Sim DW, Oh JJ, Kim KS, Choi YD, Kwon YS. Pleomorphic adenoma of the trachea. *J Bronchology Interv Pulmonol*. 2014 Jul; 21(3): 230 – 3.
- 19 Reichle G. Die Argonplasma-Koagulation zur bronchoskopischen Rekanalisation und Blutstillung. *Atemw- Lungenkrkh* 2003; Jahrgang 29: 258 – 269.
- 20 Reichle G, Freitag L, Kullmann HJ, Prenzel R, Macha HN, Farin G. Argon plasma coagulation in bronchology: a new method – alternative or complementary?. *Pneumologie* 2000; 54: 508 – 516.
- 21 Chen YK, Jakribettuu V, Springer EW, Shah RJ, Penberthy J, Nash SR. Safety and efficacy of argon plasma coagulation trimming of malpositioned and migrated biliary metal stents: a controlled study in the porcine model. *Am J Gastroenterol* 2006; 101: 2025 – 2030

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