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- 2 Eickhoff A, Repici A, Manner H, Enderle, MD. *Electrosurgical Pocket Guide for GI Interventions*; Erbe Elektromedizin GmbH
- 3 Eickhoff, A et al. *Prospective nonrandomized comparison of two modes of argon beamer (APC) tumor desobstruction: effectiveness of new pulsed APC versus forced APC*; *Endoscopy* 2007
- 4 Zenker, M. *Argon plasma Coagulation*. *GMS Krankenhhyg Interdiszip*. 2008
- 5 Taghavi SA, Soleimani SM, Hosseini-Asl SM et al. *Adrenaline injection plus argon plasma coagulation versus adrenaline injection plus hemoclips for treating high-risk bleeding peptic ulcers: a prospective, randomized trial*. *Can J Gastroenterol* 2009; 23(10): 699 – 704.
- 6 Wang HM, Hsu PI, Lo GH et al. *Comparison of hemostatic efficacy for argon plasma coagulation and distilled water injection in treating high-risk bleeding ulcers*. *J Clin Gastroenterol* 2009; 43(10): 941 – 945.
- 7 Herrera S, Bordas JM, Llach J et al. *The beneficial effects of argon plasma coagulation in the management of different types of gastric vascular ectasia lesions in patients admitted for GI hemorrhage*. *Gastrointest Endosc* 2008; 68(3): 440 – 446.

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
5 Taghavi SA, Soleimani SM, Hosseini-Asl SM et al. Adrenaline injection plus argon plasma coagulation versus adrenaline injection plus hemoclips for treating high-risk bleeding peptic ulcers: a prospective, randomized trial. *Can J Gastroenterol* 2009; 23(10): 699 – 704.

6 Wang HM, Hsu PI, Lo GH et al. Comparison of hemostatic efficacy for argon plasma coagulation and distilled water injection in treating high-risk bleeding ulcers. *J Clin Gastroenterol* 2009; 43(10): 941 – 945.

7 Herrera S, Bordas JM, Llach J et al. The beneficial effects of argon plasma coagulation in the management of different types of gastric vascular ectasia lesions in patients admitted for GI hemorrhage. *Gastrointest Endosc* 2008; 68(3): 440 – 446.

8 Manner H, Enderle MD, Pech O et al. Second-generation argon plasma coagulation: two-center experience with 600 patients. *J Gastroenterol Hepatol* 2008; 23(6): 872 – 878.

9 Katsinelos P, Gkagkalis S, Paroutoglou G, Chatzimavroudis G, Fasoulas K, Zavos C, Varitimadakis K, Lazaraki G, Kotronis G, Kountouras J. A prospective comparative study of blended and pure coagulation current in endoscopic mucosal resection of large sessile colorectal polyps. *Surg Laparosc Endosc Percutan Tech*. 2014 Jun; 24(3): 226 – 31.

10 Manner H, Rabenstein T, Pech O, Braun K, May A, Pohl J, Angelika Behrens A, Vieth M, Ell C. Ablation of residual Barrett's epithelium after endoscopic resection: a randomized long-term follow-up study of argon plasma coagulation vs. surveillance (APE study). *Endoscopy* 2014; 46(01): 6–12

15 Wang H, Tao M, Zhang N, Luo L, Li D, Zou H, Zhou Y, Liang S. Bronchoscopic interventions combined with percutaneous modalities for the treatment of thyroid cancers with airway invasion. *Eur Arch Otorhinolaryngol*. 2015 Feb; 272(2): 445 – 51.

16 Wang JW, Huang M, Zha WJ, Zhou LF, Qi X, Wang H. Flexible bronchoscopic intervention for endobronchial hamartoma. *Zhonghua Jie He Hu Xi Za Zhi*. 2013 Dec; 36(12): 963 – 7. Chinese


17 Seaman JC, Musani AI. Endobronchial ablative therapies. *Clin Chest Med*. 2013 Sep; 34(3): 417 – 25. doi: 10.1016/j.ccm.2013.04.006.

18 Sim DW, Oh JJ, Kim KS, Choi YD, Kwon YS. Pleomorphic adenoma of the trachea. *J Bronchology Interv Pulmonol*. 2014 Jul; 21(3): 230 – 3.

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21 Chen YK, Jakribettuu V, Springer EW, Shah RJ, Penberthy J, Nash SR. Safety and efficacy of argon plasma coagulation trimming of malpositioned and migrated biliary metal stents: a controlled study in the porcine model. *Am J Gastroenterol* 2006; 101: 2025 – 2030

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Important information

We have prepared this document with care. Nonetheless, we cannot completely rule out errors in this document.

The information, recommendations and other data („Information“) contained in this document reflect our state of knowledge and the state of science and technology at the time of preparing the document. The information is of a general nature, non-binding and serves solely for general information purposes and does not represent instructions for use or notes on application.

The information and recommendations contained in this document do not constitute any legal obligations by Erbe Elektromedizin GmbH as well as their associated companies („Erbe“) or any other claims against Erbe. The information does not represent a guarantee or other quality statement, these require the express contractual arrangement with Erbe in individual cases. Erbe shall not be liable for any type of damage resulting from following information given in this document, regardless of the legal reason for liability.

Every user of an Erbe product is responsible for checking the respective Erbe product for its properties as well as the suitability for the intended type of application or intended purpose in advance. The suitable type of application of the respective Erbe product is given by the user manual and the notes on use for the corresponding Erbe product. The user is obliged to check whether the existing user manual and the notes on use correspond with the status for the specific Erbe product. The devices may only be used according to the user manual and the notes on use.

The information on setting values, application sites, duration of application and the use of the respective Erbe product is based on the clinical experience of physicians independent from Erbe. They represent guidelines which need to be checked by the user for their suitability for the actual planned application. Depending on the circumstances of an actual application case, it may be necessary to deviate from the information provided. The user has to check this on his/her own responsibility in each case when using an Erbe product. We wish to point out that science and technology is constantly subject to new developments arising from research and clinical experience. For this reason it may be necessary for the user to deviate from the information provided in this document.

This document contains information about Erbe products which may possibly not be approved in a specific country. The user of the respective Erbe product is obliged to inform him/herself whether the Erbe product he/she is using is legally approved in his/her country and/or if legal requirements or restrictions for use possibly exist and to which extent.

This document is not intended for users in the USA.

Erbe Elektromedizin GmbH
Waldhoernlestrasse 17
72072 Tuebingen
Alemania

Tel +49 7071 755-0
Fax +49 7071 755-179
info@erbe-med.com
erbe-med.com