

# Precise application in complex thin-walled anatomy

*Benefit from an all-in-one solution uniting cutting-edge technology with new flexibility.*

*Inject, cut, coagulate – no need for instrument exchange – for efficient workflow and ease of use.*

**TAILOR YOUR OWN TECHNIQUE.**



**HYBRIDknife® flex  
I-Type**

Elevating  
3rd space endoscopy

**erbe**  
power your performance.

Erbe Elektromedizin GmbH · Waldhoernlestrasse 17 · 72072 Tuebingen · Germany · +49 7071 755-0 · [info@erbe-med.com](mailto:info@erbe-med.com)

JOIN US ON     





# Elevating 3rd space endoscopy

**Thin electrode** for precise applications with limited thermal damage

**Coaxial nozzle** for high-pressure injection

**Flexible probe** for smooth passage even in retroflexion

**Small diameter** for adequate suction through the working channel

**Locking mechanism** for constant electrode length

## Inject needle-free

Expand the submucosal space without instrument exchange using high-pressure, needle-free, tissue-selective hydrodissection.

## Cut and Coagulate precisely

Achieve high electrosurgical precision with limited thermal damage using the thin electrode, and apply various coagulation techniques with VIO® modes.

## Flex with stability

Expect smooth passage and stable electrode length, even in difficult angled or retroflexed endoscopic positions.

Marketing claim list 0188716, V005

erbe

VIO® 3

APC3

ERBEJET® 2

**THE HYBRIDKNIFE® FLEX  
UNITES ELECTROSURGERY WITH  
HYDRODISSECTION**

Powered by VIO® 3



Tubing length	2300 mm
Outer diameter (distal end)	2.6 mm
Outer diameter (tubing)	2.1 mm
Electrode length	1.5 mm / 2.0 mm
Type	I
Retro-flexibility	More than 210°
I-Type 1.5 mm	No. 20150-112
I-Type 2.0 mm	No. 20150-113

FIND YOUR PRODUCT